










MIS SEVEN

“Crown Down” drilling chart

Implant	Pilot drill L= 5 mm	Final drill Final Length
 3.30 NP	3.30 mm	2.50 mm
 3.75 SP	3.75 mm	3.30 mm
 4.20 SP	4.20 mm	3.75 mm
 5.00 WP	5.00 mm	4.50 mm
 6.00 WP	6.00 mm	5.50 mm

MIS C1

“Crown Down” drilling chart

Implant	Pilot drill L= 5 mm	Final drill Final Length
 3.30 NP	3.30 mm	2.50 mm
 3.75 SP	3.75 mm	3.30 mm
 4.20 SP	4.20 mm	3.75 mm
 5.00 WP	5.00 mm	4.50 mm