

# “Crown Down” drilling chart

Neo Biotech®

| Implant   | Pilot drill<br>L= 3-5 mm | Final drill<br>Final Length |
|---|--------------------------|-----------------------------|
|  3.5   | 3.75mm                   | 3.00 mm                     |
|  4.0  | 4.2 mm                   | 3.75 mm                     |
|  4.5 | 4.50 mm                  | 4.00mm                      |
|  5.0 | 5.00 mm                  | 4.5 mm                      |
|  5.5 | 5.50 mm                  | 5.00mm                      |
|  6.0 | 6.00 mm                  | 5.50 mm                     |