








Straumann® BLX

“Crown Down” drilling chart

Implant	Pilot drill L= 5 mm	Final drill Final Length
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px;">Ø 3.5</div>  </div>	Ø 3.5	Ø 3.0
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px;">Ø 3.75</div>  </div>	Ø 3.75	Ø 3.3
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px;">Ø 4.0</div>  </div>	Ø 4.0	Ø 3.5
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px;">Ø 4.5</div>  </div>	Ø 4.5 Ø 3.5 D3-D4 BONE	Ø 4.0 Ø 3.5 D3-D4 BONE
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px;">Ø 5.0</div>  </div>	Ø 5.0	Ø 4.5
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px;">Ø 5.5</div>  </div>	Ø 5.5 Ø 4.5 D3-D4 BONE	Ø 5.0 Ø 4.5 D3-D4 BONE
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 5px;">Ø 6.5</div>  </div>	Ø 6.0 Ø 5.5 D3-D4 BONE	Ø 5.5 Ø 5.5 D3-D4 BONE